


PHILIPS



Resuscitation Supplies and Accessories Catalog

2007 Edition

<http://shop.medical.philips.com>

HEARTSTART
DEFIBRILLATORS

Defibrillator Supplies/Accessories

External Pads and Cables	1
Pad Adapters	3
Adapter Cable Test Loads	3
External Paddle Sets	4
Replacement Electrodes	4
Internal Paddles	5
Sync Cables	6

ECG Cables and Leads

Patient Cables and Leads	7
--------------------------------	---

Miscellaneous Supplies and Accessories

ECG Monitoring Electrodes	11
Recording Paper	11
Defibrillator Batteries/Cases/Accessories	12
Data Management Solutions	13

Non-Invasive Blood Pressure

Reusable Non-Invasive Blood Pressure Cuffs	15
Disposable Non-Invasive Blood Pressure Cuffs	17
Interconnect Tubing	18

SpO₂ Monitoring

Reusable SpO ₂ Sensors	19
Disposable SpO ₂ Sensor	20

EtCO₂ Monitoring

Microstream CO ₂	21
-----------------------------------	----

Invasive Pressure and Continuous Temperature

Invasive Pressure	23
Continuous Temperature	24

AED Supplies and Accessories

AED Pads and Accessories	25
Data Management Solutions	28
AED Cases and Wall Mount Solutions	30
AED Training Materials	32

Using This Guide

Using This Guide

This guide lists most medical supplies for *Resuscitation Monitoring*. All information in this guide is subject to change without notice. For a comprehensive list of medical supplies across all departments, refer to the following sources:

- *Medical Supplies Reference Guide* (Contact your local Sales Rep for the latest edition.)
- The Medical Supplies product catalog, which is updated daily: <http://shop.medical.philips.com>

Philips Medical Systems Medical Supplies Guarantee

Philips Medical Systems is committed to your satisfaction. That's why we guarantee virtually all of our products against defects in materials and workmanship.

Ordering Information

To order the medical supplies in this guide or to obtain more information or additional literature, contact your local Medical Supplies representative or visit our website: <http://shop.medical.philips.com>.

Online Ordering (U.S. Only)

Online ordering is available to U.S. customers at <http://shop.medical.philips.com>. You can search for supplies by application category, Philips instrument, or keyword.

1. Log in using a secure password.
2. Enter the product number and quantity in the shopping cart.
3. Enter your PO or credit card number at check-out.
4. Your order is immediately transferred into our fulfillment system, where you can track it anytime.

Latex Cautions

Philips Medical Systems currently distributes over 500 products, of which less than 9% contain natural rubber latex or dry natural rubber that could contact humans. Because these substances may cause allergic reactions, products in this guide that contain them are identified with the following symbol:




Caution: Products that are identified with this symbol contain dry natural rubber or natural rubber latex that can come into contact with humans and may cause allergic reactions.

To the best of our knowledge, products not identified with this symbol are not known to contain natural rubber latex or dry natural rubber.

The latex symbol indicates **only** which **medical supplies** contain natural rubber latex or dry natural rubber that can come into contact with humans and may cause allergic reactions. If you want to know whether a Philips **medical device** contains natural rubber latex or dry natural rubber, refer to the product labeling.

The FDA requires that products be labeled with human contacting latex content information when appropriate. Always refer to product labeling for accurate information on human contacting dry natural rubber or natural rubber latex.



CE Marking

Products labeled with this symbol  conforms to the European Union (EU) directives as mandated by EU member state legislation and are available for sale in Europe to countries that accept CE marking with restrictions as may be noted. Medical devices labeled with this marking means the product is in conformity with the European Directive (MDD): 93/42/EEC Council Directive of 14 June 1993 on the approximation of the laws of the Member States concerning medical devices.

In order for a product to be able to be labeled with the CE marking, product quality and manufacturing quality systems must be documented to comply with certain standards. Compliance with the European directives and regulations ensures a standard level of performance for device function, safety, and health protection for users, patients, and others.

Non-medical products that comply with other harmonized directives may also bear the CE mark.

Countries that accept CE marked products include Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Liechtenstein, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Russia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and the United Kingdom. A CE mark does not necessarily indicate that a product is marketed or distributed in all of the countries listed above.

Each Philips medical product that is labeled with the Philips Medical Systems CE marking will be identified with  within this guide. Philips'  marking can be indicative of any of the the following:



Other manufacturers' CE marking symbols are indicated below:






Defibrillator Supplies/Accessories

External Pads and Cables

- For non-invasive pacing, hands-free defibrillation, synchronized cardioversion, and ECG monitoring

Disposable External Multifunction Electrode Pads

Product No.	Description	Shelf Life	Sets/Case	For Philips Instrument
M3713A CE	HeartStart Adult Plus Pads with plug-style connector	24 mths from date of mfg.	10	E, S, EM (Heartstream ForeRunner AED), 940010XX, 940020XX, 940030XX (Heartstart FR AED), M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 AED) M3500B (Heartstream/HeartStart XLT) M3535A/M3536A/M3536j (HeartStart MRx/MRxE) M4735A (Heartstream/HeartStart XL) M5500B (Heartstart 4000)
M3716A CE	HeartStart Adult Radiolucent Pads with plug style connector			
M3718A CE	HeartStart Adult Radiotransparent/Reduced Skin Irritation Pads with plug style connector			
M3717A CE	HeartStart Pediatric Plus Pads with plug style connector			
M3719A CE	HeartStart Pediatric Radiotransparent/Reduced Skin Irritation Pads with plug style connector	5	5	M3500B (Heartstream/HeartStart XLT) M3535A/M3536A/M3536j (HeartStart MRx/MRxE) M4735A (Heartstream/HeartStart XL) M5500B (Heartstart 4000)
M3501A AAMI	Philips Adult Pads with barrel style connector			
M3502A CE IEC				
M3503A CE IEC				
M3504A AAMI				



M3716A



M3717A



M3718A

External Multifunction Cables

Product No.	Description		For Philips Instrument
M3507A CE	Hands-free Pad cable	For use with Philips pads with barrel style connector. Can also be used with plug style connector pads with 05-10200 plug-to-barrel style connector adapter. Length: 2.1m (7 ft.)	M3500B (Heartstream/HeartStart XLT) M3535A/M3536A (HeartStart MRx)
M3508A CE	HeartStart Hands-free Pad cable	For use with Philips pads with plug style connector. Length: 2.1m (7 ft.)	M4735A (Heartstream/HeartStart XL) M5500B (Heartstart 4000)
M4748A CE	HeartStart extension cable	Length = 10cm (4 in)	M4735A (Heartstream/HeartStart XL) M5500B (Heartstart 4000)



M3507A



M3508A

Pad Adapters


HeartStart Pad Adapters

Product No.	Description	Use with Instrument Models
05-10000 CE	Removable adapter	Medtronic Physio-Control LifePak 500, 9, 10c, 11, 12, 20
05-10100 CE	Removable adapter	Zoll 1200, 1400, 1600, and M-Series
05-10200 CE	Removable or permanent adapter	Fits Philips hands-free cable M3507A for connection to: M4735A (Heartstream/HeartStart XL); M5500B (Heartstart 4000) M3535A/M3536A (HeartStart MRx)


Adapter Cable Test Loads

Test Loads


Product No.	Description
M1781A CE	50-ohm test load for M3507A Multifunction cables.
M3725A CE	HeartStart 50-ohm test load for the M3508A HeartStart Hands-Free cable. Also works with M3507A with 05-10200 plug-to-barrel style connector adapter.



05-10000



05-10100



05-10200

External Paddle Sets

External Anterior/Anterior Paddle Sets

Product No.	For Philips Instrument	Sterilization
M3542A With PCI CE	M3535A/M3536A (HeartStart MRx)	Not designed to be sterilized. Refer to operating manual instructions for cleaning exterior surfaces.
M3543A With PCI CE (water resistant)	M3535A/M3536A (HeartStart MRx)	Not designed to be sterilized. Refer to operating manual instructions for cleaning exterior surfaces.
M4745A Without PCI CE	M4735A (Heartstream/HeartStart XL) M3535A/M3536A (HeartStart MRx)	Clean and sterilize according to operating manual instructions
M4746A With PCI CE	M4735A (Heartstream/HeartStart XL)	Not designed to be sterilized. Refer to operating manual instructions for cleaning exterior surfaces.

Replacement Electrodes

Replacement Electrodes for External Paddle Sets

Product No.	Description	Philips External Paddles
M1789A CE	Rectangular Paddle Electrode Replacement (clip-on)	M4746A
M4759A	Rectangular Paddle Electrode Replacement (clip-on)	M3542A



M3542A


Internal Paddles

Switched Reusable Internal Paddles (Discharge Switch on Handle)

Product No.	Diameter	Sterilization	For Philips Instrument
M4741A CE	7.5cm (3in) Adult	Clean and sterilize according to operating manual instructions	M4735A (Heartstream/HeartStart XL) M3535A/M3536A (HeartStart MRx)
M4742A CE	6cm (2.4in) Adult		
M4743A CE	4.5cm (1.8in) Pediatric		
M4744A CE	2.8cm (1.1in) Infant		

Switchless Reusable Internal Paddles (Discharge Switch on Defibrillator)

Product No.	Diameter	Sterilization	Use Adapter Cable	For Philips Instrument
M1741A CE	7.5cm (3in) Adult	Clean and sterilize according to operating manual instructions	For M4735A and M3535A/M3536A defibrillators use M4740A adapter cable	M4735A (Heartstream/HeartStart XL) M3535A/M3536A (HeartStart MRx)
M1742A CE	6cm (2.4in) Adult			
M1743A CE	4.5cm (1.8in) Pediatric			
M1744A CE	2.8cm (1.1in) Infant			



M4741A-M4744A



M1741A-M1744A

Internal Paddle Adapters

Product No.	Description	Length	Use With Internal Paddle Set	For Philips Instrument
M4740A CE	HeartStart Switchless Internal Paddle Adapter	20cm (8 in)	Switchless Reusable Internal Paddles M174xA series	M4735A (Heartstream/HeartStart XL) M3535A/M3536A (HeartStart MRx)
M4748A CE	HeartStart Extension Cable	10cm (4 in)	Switched Reusable Internal Paddles M474xA series	M4735A (Heartstream/HeartStart XL)



M4740A

Sync Cables

Sync Cables

Product No.	Connector	For Philips Defibrillators	Length
M1783A CE	12-pin	M4735A (Heartstream/HeartStart XL) M3535A/M3536A/M3536J (HeartStart MRx/MRxE)	2.5m (8ft)
M5526A CE	12-pin	M3535A/M3536A/M3536J (HeartStart MRx/MRxE) M4735A (Heartstream/HeartStart XL)	7.8m (25ft)



ECG Cables and Leads

Patient Cables and Leads

3-Lead Monitoring

Number of Pins in Connector to Defib	Color Code	Order Trunk Cable	with Lead Set	For Philips Instrument
8-pin	AAMI	M1733A	N/A (one piece cable and leads)	M3500B (Heartstream/HeartStart XLT) M5500B (Heartstart 4000)
	IEC (Europe)	M1735A CE	N/A (one piece cable and leads)	
	AAMI	M1580A	M1603A (Grabber) or M1605A (Snap)	M3500B (Heartstream/HeartStart XLT) M5500B (Heartstart 4000)
	IEC (Europe)	M1590A CE	M1613A (Grabber) or M1615A (Snap)	
12-pin	AAMI	M1500A	M1603A (Grabber) or M1605A (Snap)	M3535A/M3536A/M3536j (HeartStart MRx/MRxE) M4735A (Heartstream/HeartStart XL)
	IEC (Europe)	M1510A CE	M1613A (Grabber) or M1615A (Snap)	
	AAMI	M1669A	M1671A (Grabber) or M1673A (Snap)	M3535A/M3536A/M3536j (HeartStart MRx/MRxE) M4735A (Heartstream/HeartStart XL)
	IEC (Europe)	CE	M1672A (Grabber) or M1674A (Snap)	




M1500A



M1669A



M1671A



M1672A

5-Lead Monitoring

Number of Pins in Connector to Defib	Color Code	Order Trunk Cable	with Lead Set	For Philips instrument
8-pin	AAMI	M1734A	N/A (one piece cable and leads)	M3500B (Heartstream/HeartStart XLT) M5500B (Heartstart 4000)
	IEC (Europe)	M1736A CE	N/A (one piece cable and leads)	
	AAMI	M1600A	M1623A (Grabber) or M1625A (Snap)	M3500B (Heartstream/HeartStart XLT) M5500B (Heartstart 4000)
	IEC (Europe)	M1610A CE	M1633A (Grabber) or M1635A (Snap)	
12-pin	AAMI	M1520A	M1623A (Grabber) or M1625A (Snap)	M3535A/M3536A/M3536j (HeartStart MRx/MRxE) M4735A (Heartstream/HeartStart XL)
	IEC (Europe)	M1530A CE	M1633A (Grabber) or M1635A (Snap)	
	AAMI	M1668A	M1968A (Grabber) or M1644A (Snap)	M3535A/M3536A/M3536j (HeartStart MRx/MRxE) M4735A (Heartstream/HeartStart XL)
	IEC (Europe)	CE	M1971A (Grabber) or M1645A (Snap)	




M1520A



M1668A



M1971A




M1968A


10-Electrode Cable Sets - 5 Wire

- For 10-Lead applications, choose one 'chest' and one 'limb' lead set

Number of Pins in Connector to Defib	Color Code	Order Trunk Cable	Lead Set	For Philips Instrument
12-pin	AAMI	M1663A	M1968A ICU (Grabber) or M1644A (Snap) limb	M3535A/M3536A/M3536J (HeartStart MRx/MRx E)
			M1976A ICU (Grabber) or M1602A (Snap)- chest	
	IEC (Europe)	M1663A CE	M1971A ICU (Grabber) or M1645A (Snap) - limb	
			M1978A ICU (Grabber) or M1604A (Snap)- chest	




M1971A



M1968A



M1978A



M1976A



M1663A

10-Electrode Cable Sets - 3 and 7 Wire

Number of Pins in Connector to Defib	Color Code	Order Trunk Cable	with Lead Set	For Philips Instrument
12-pin	AAMI	M3525A	M3526A 3-wire (Snap)	M3535A/M3536A/M3536J (HeartStart MRx/MRxE)
			M3527A 7-wire (Snap) (for 12-lead use with M3526A)	
	IEC (Europe)	M3525A CE	M3528A 3-wire (Snap)	
			M3529A 7-wire (Snap) (for 12-lead use with M3528A)	
Combiner Plug		M5530A CE	For 3-wire lead set M3526A and M3528A	



M3525A/M3526A/M3527A



Miscellaneous Supplies and Accessories

ECG Monitoring Electrodes

- Disposable radiolucent electrodes for use where x-rays will be used

Disposable Radiolucent Electrodes

Product Number	Description	Dimensions	Packaging	Shelf life:	Philips Instrument
M2202A CE See photo	Silver/silver-chloride (Ag/AgCl) sensor; non-metallic snap, pre-gelled.	2" (50mm) diameter	5/pack 1 case= 60 packs= 300 electrodes	At least 6 months	Heartstream/HeartStart XLT Heartstream/HeartStart XL Heartstart 4000 HeartStart MRx/MRxE

Recording Paper

- Designed to minimize maintenance and repair costs while providing optimum performance for your recorder

Chemical/Thermal Recording Paper

Product No.	Format	Grid Width	Critical Margin	Qty/Case	Philips Instrument
40457C CE	1-channel roll	40mm (1.6in)	6mm (0.25in)	10 rolls	M3500B (Heartstream/HeartStart XLT) M3535A/M3536A/M3536J (HeartStart MRx/MRxE)
40457D CE				80 rolls	
989803138171 CE		65mm (2.5in)	10mm (0.4in)	10 rolls	M3535A/M3536A/M3536J Opt. B04 75mm printer (HeartStart MRx/MRxE)
989803138181 CE				80 rolls	



M2202A



40457C

Defibrillator Batteries/Cases/Accessories

Defibrillator Batteries

Product No.	For Philips Instrument	Description	Shelf Life
M3506A CE	M3500B (Heartstream/HeartStart XLT) M4735A (Heartstream/HeartStart XL) M5500B (Heartstart 4000)	Battery charger adapter	N/A
M3516A CE	M3500B (Heartstream/HeartStart XLT) M4735A (Heartstream/HeartStart XL) M5500B (Heartstart 4000)	Sealed 12-volt lead acid (SLA) battery	At least 6 months
M3538A CE	M3535A/M3536A/M3536J (HeartStart MRx/MRxE)	14.8-volt Lithium Ion battery	N/A



M3516A



M3538A

Defibrillator Cases and Accessories

Product No.	Description	For Philips Instrument
M3509A	Carrying case	M3500B (Heartstream/HeartStart XLT)
M5509A	Carrying case	M5500B (HeartStart 4000)
M4751A	Accessory pouch (red)	M4735A (HeartStart XL)

Data Management Solutions

HeartStart Software and Data Card

Product No.	Description	For Philips Instrument
861276 CE	HeartStart Event Review Pro Software Full featured software that allows you to review and create reports on sudden cardiac arrest event data recorded by Philips defibrillators.	M3854A, M3524A (FR2 series data card and data card reader) M3500B (Heartstream/HeartStart XLT) M4735A (Heartstream/HeartStart XL) M3535A, M3536A, M3536J (HeartStart MRx/MRxE) M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 series AED) E, S, EM, (Heartstream ForeRunner AED) 940010XX, 940020XX, 940030XX (Heartstart FR AED) M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)
989803142521 CE	HeartStart 12-lead Transfer Station. Software product received a 12-lead ECG transmission from an appropriately equipped HeartStart MRx defibrillator, archives it and reliably forwards it to the appropriate hospital destination.	M3535A, M3536A/M3536J (HeartStart MRx/MRxE)
989803147711 CE	HeartStart data card (one card per order)	M4735A (Heartstream/HeartStart XL)
989803147721 CE		M3500B (Heartstream/HeartStart XLT) M5500B (HeartStart 4000)

HeartStart MRx Accessories

Product No.	Description	Philips Instrument
M3537A CE	Bed Rail Hook	M3535A/M3536A/M3536J (HeartStart MRx/MRxE Monitor/Defibrillator)
M3541A CE	Carrying case (3 accessory pouches and shoulder strap)	
M3544A CE	Data Card Tray	
M3545A CE	Data Card and Tray	
M5521A CE	Color handle (Green)	
M5522A CE	Color handle (Blue)	
M5523A CE	Color handle (Yellow)	
M5524A CE	Color handle (Rose)	
M5525A CE	Color handle (Grey)	
M5528A	Vehicle Wall Mount	
M5529A	DC Power Module	



M5521A-M5525A



M3545A



M3541A



M5528A



M5529A



Non-Invasive Blood Pressure

Reusable Non-Invasive Blood Pressure Cuffs

Antimicrobial Coated Multi-Patient NIBP Cuffs

- One year warranty

Single Hose Product No.	Cuff Size/ Color	Limb Circumference	Philips Instrument	Qty	Weight
M4552A CE	Infant (Orange)	Bladder width=2.1" (5.4cm) Bladder length=4.8" (12.3cm) Limb Circumference= 9-15cm	M3535A/M3536A (HeartStart MRx)	1 cuff per bag	255.1g (9oz)
M4553A CE	Pediatric (Green)	Bladder width=3.1" (8cm) Bladder length=7.1" (18.1cm) Limb Circumference = 15-21.5cm			425.2g (15oz)
M4554A CE	Small Adult (Light Blue)	Bladder width= 4.2" (10.6cm) Bladder length= 9.4" (23.9cm) Limb Circumference = 20.5-28.5cm			581.17g (1.3lb)
M4555A CE	Adult (Blue)	Bladder width=5.3" (13.5cm) Bladder length=12.1" (30.7cm) Limb Circumference = 27.5-36.5cm			734.25g (1.6lb)
M4557A CE	Large Adult (Dark Red)	Bladder width=6.7" (17cm) Bladder length=15.2" (38.6cm) Limb Circumference = 35.5-46cm			1.07kg (2.4lb)
M4559A CE	Thigh (Grey)	Bladder width=8.3"(21cm) Bladder length=18.7" (47.5cm) Limb Circumference = 45-56.5cm			1.28kg (2.8lb)



Antimicrobial Multi-Patient NIBP Cuffs (Color Coded Sizes)

Reusable Comfort Cuffs

- One year warranty

Product No.	Description	Limb Circumference	Philips Instrument	Qty	Weight
M1572A CE	Pediatric size cuff	Bladder width=3.1" (8cm) Bladder length=8.5" (21.5cm) Limb circumference = 14cm to 21.5cm	M3535A/M3536A (HeartStart MRx)	1 cuff	40g (1.4oz)
M1573A CE	Small adult size cuff	Bladder width=4.1" (10.5cm) Bladder length=9.4" (24cm) Limb circumference= 20.5cm to 28cm			57g (2oz)
M1574A CE	Adult size cuff	Bladder width=5.1" (13cm) Bladder length=11" (28cm) Limb circumference = 27cm to 35cm			74g (2.6oz)
M1575A CE	Large adult size cuff	Bladder width=6.3" (16cm) Bladder length=13.8" (35cm) Limb circumference = 34cm to 43cm			106g (3.7oz)

Reusable Traditional Cuffs


- 3-Year Warranty

Product No.	Description	Limb Circumference	Philips Instrument	Qty/ Pack	Weight
40401A CE	Infant size cuff	Bladder width=2.4" (6cm) Bladder length=4.7" (12cm) Limb circumference= 10cm to 19cm	M3535A/M3536A (HeartStart MRx)	1 cuff	57g (2oz)
40401B CE	Pediatric size cuff	Bladder width=3.5" (8.9cm) Bladder length=7.1" (18cm) Limb circumference= 18cm to 26cm			85g (3oz)
40401C CE	Adult size cuff	Bladder width= 4.8"(12.1cm) Bladder length=9.1" (23cm) Limb circumference= 25cm to 35cm			142g (5oz)
40401D CE	Large arm size cuff	Bladder width=6"(15.2cm) Bladder length=13"(33cm) Limb circumference= 33cm to 47cm			197g (7oz)
40401E CE	Thigh size cuff	Bladder width=7"(17.8cm) Bladder length=14.2"(36cm) Limb circumference= 46cm to 66cm			255g (9oz)

Disposable Non-Invasive Blood Pressure Cuffs

Single-Patient Disposable NIBP Cuffs

Single Hose Product Number	Size (Print Color)	Description	Use with Philips Instrument	Qty	Weight
M4572A CE	Infant (Orange)	Bladder width=2.1" (5.4cm) Bladder length=4.8" (12.3cm) Limb Circumference= 9-15cm	M3535A/M3536A (HeartStart MRx)	10 per pack	Pkg wt=9oz (255.1g)
M4573A CE	Pediatric (Green)	Bladder width=3.1" (8cm) Bladder length=7.1" (18.1cm) Limb Circumference = 15-21.5cm			Pkg wt=15oz (425.2g)
M4574A CE	Small Adult (Light Blue)	Bladder width= 4.2" (10.6cm) Bladder length= 9.4" (23.9cm) Limb Circumference = 20.5-28.5cm			Pkg wt=19oz (581.17g)
M4575A CE	Adult (Blue)	Bladder width=5.3" (13.5cm) Bladder length=12.1" (30.7cm) Limb Circumference = 27.5-36.5cm			Pkg wt=22oz (734.25g)
M4576A CE	Adult Extra Long (Blue)				Pkg wt=2.4lb (1.1kg)
M4577A CE	Large Adult (Dark Red)	Bladder width=6.7"(17cm) Bladder length=15.2" (38.6cm) Limb Circumference = 35.5-46cm			Pkg wt=2.4lb (1.07kg)
M4578A CE	Large Adult Extra Long (Dark Red)				Pkg wt=2.7lb (1.2kg)
M4579A CE	Thigh (Brown)	Bladder width=8.3"(21cm) Bladder length=18.7" (47.5cm) Limb Circumference = 45-56.5cm			Pkg wt=2.8lb (1.28kg)



Soft Single-Patient Disposable NIBP Cuffs

Interconnect Tubing

NIBP Interconnect Tubing

- For use with adult/pediatric NIBP cuffs only

Product No.	Description	Length	Qty/Case	Weight
M1598B CE	Connects Philips's reusable and disposable adult/pediatric NIBP cuffs to monitor.	1.5m (4.9ft)	1 bag (1 tubing)	85g (3oz)
M1599B CE	Connects Philips's reusable and disposable adult/pediatric NIBP cuffs to monitor. (Not available in Japan)	3m (9.8ft)		142g (5oz)



M1599B



SpO₂ Monitoring

Reusable SpO₂ Sensors


Reusable SpO₂ Sensors and Cables

- Reusable sensors cut costs with greater ease of use
- M119X series sensors are sanitized by immersing in cleaning solution


Product No.	Description	Philips Instrument
M1191A CE	Adult finger SpO ₂ sensor Patient size >50kg (>110 lb)	M3500B (HeartStart XLT)
M1192A CE	Pediatric finger SpO ₂ sensor Patient size 15–50kg (33–110 lb)	M3535A/M3536A/M3536J/ (HeartStart MRx/MRxE)
M1194A	Adult ear clip SpO ₂ sensor	M4735A (Heartstream/HeartStart XL, XLT)
M1195A	Infant	M5500B (HeartStart 4000)
Cable		
M1941A CE	Extension cable for Philips reusable SpO ₂ sensors (1 cable per pack). Use with M3500B, M3535A/M3536A/M3536J, M4735A, M5500B defibrillators. Length: 6.5 ft. (2m)	



M1191A



M1192A




M1941A

Disposable SpO₂ Sensor

- Pre-aligned LED
- One size fits Adults and Pediatrics
- Designed to perform for a single patient's average hospital stay

Product No.	Description	Philips Instrument	Adapter Cable
M1131A CE	Adult/Pediatric disposable SpO ₂ sensor Patient size: >20kg (>44lb.)	M3500B (HeartStart XLT) M5500B (HeartStart 4000) M3535A/M3536A/M3536J (HeartStart MRx/MRxE) M4735A (Heartstream/HeartStart XL)	Use M1943A (1.1m)



M1131A

Adapter Cable

Product No.	Description
M1943A CE	Adapter Cable (8 pin to 9 pin D-sub connector) Length: 1.1m



M1943A




EtCO₂ Monitoring

Microstream CO₂

Use with M3535A/M3536A Opt. A03 and A04 (HeartStart MRx)

Ventilation	Environment	Patient Weight	Accessories	Qty/ Pack	Order Number
Intubated	Non-humidified	>2kg	FilterLine Set Adult/Pediatric	25 per case	M1920A CE
		<2kg	FilterLine H Set Infant/Neonatal		M1923A CE
	Humidified	>2kg	FilterLine H Set Adult/Pediatric		M1921A CE
		<2kg	FilterLine H Set Infant/Neonatal		M1923A CE
Non-Intubated, Single Purpose	Oral-Nasal CO ₂ Non-Humidified	>20kg	Smart CapnoLine Plus Adult		M2526A CE
		10-20kg	Smart CapnoLine Pediatric		M2524A CE
Non-Intubated, Dual Purpose	Oral-Nasal CO ₂ + O ₂ Non-Humidified	>20 kg	Smart CapnoLine Plus O ₂ Adult		M2522A CE
		10-20kg	Smart CapnoLine O ₂ Pediatric		M2520A CE





Invasive Pressure and Continuous Temperature

Invasive Pressure

Memscap Reusable Pressure Transducer

Product No.	Description	Qty.	Philips Instrument
CPJ840J6 CE	Reusable pressure transducer	1	M3535A/M3536A Opt. A06
CPJ84022 CE	Sterile disposable pressure domes for use with CPJ840J6	50	
CPJ84046 CE	Transducer holder for CPJ840J6	4	

Disposable Invasive Pressure Kits (available in Europe and Asia only)


Product No.	Description	Qty.	Philips Instrument
M1567A CE	Single channel disposable blood pressure kit	1 box = 20 kits	M3535A/M3536A Opt. A06
M1568A CE	Dual line blood pressure kit for measuring CVP, ABP and other pressure measurements.		
M1634A CE	Adapter cable for M1567A and M1568A kits.	1 cable	



CPJ840J6 and CPJ84022 dome



M1567A




M1568A


Continuous Temperature

Sterile Disposable Temperature Probes

Product No.	Description	Qty.	Philips Instrument
21090A CE	Disposable esophageal/rectal temperature probe (12FR)	20/box	M3535A/M3536A Opt. A06
M1837A CE	Disposable esophageal/rectal temperature probe (9FR)		
21091A CE	Disposable skin surface temperature probes (28" long)		
21093A CE	Disposable esophageal stethoscope temperature probe (12FR)		
21094A CE	Disposable esophageal stethoscope temperature probe (18FR)		
21095A CE	Disposable esophageal stethoscope temperature probe (24FR)	10/box	
21096A CE	Disposable Foley catheter temperature probe (16FR)		
21097A CE	Disposable Foley catheter temperature probe (18FR)		
M2255A CE	Disposable Foley catheter temperature probe (14FR)		
Temperature Probe Extension Cables for Disposable Probes			
21082A CE	2-pin plug extension cable, accepts mini (1/8") phone plug. Terminates in 2-prong plug.	Length: 10'	M3535A/M3536A Opt. A06
21082B CE		Length: 5'	



21091A



21096A

Reusable Temperature Probes

Product No.	Description	Qty.	Philips Instrument
21075A CE	Esophageal/rectal, adult, 12 FR (length: 10')	1	M3535A/M3536A Opt. A06
21076A CE	Esophageal/rectal, pediatric, 10 FR (length: 10')	1	
21078A CE	Skin surface (length: 10')	1	



21075A



21078A




AED Supplies and Accessories

AED Pads and Accessories

External AED Pads

Product No.	Description	Quantity	Philips Instrument
DP2 CE	Self-adhesive adult electrode pads with integrated cable and connector for patients weighing 55 lbs. (25kg) or more or 8 years and above.	2 sets	E, S, EM (Heartstream ForeRunner AED) 940010XX, 940020XX, 940030XX (Heartstart FR AED)
DP6 CE	Self-adhesive adult electrode pads with integrated cable and connector for patients weighing 55 lbs. (25kg) or more or 8 years and above.	6 sets	M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 Series AED)
M3870A CE	Self-adhesive FR2 series infant/child electrode pads. Reduced energy (50 Joules) providing more appropriate dosage for infants and children under 8 years old or 55 lbs (25kg).	1 set	For use with the Heartstream and HeartStart FR2 Series AED models only M3860A, M3861A, M3840A, M3841A
M5071A CE	Adult SMART Pads Cartridge for patients weighing 55 lbs (25kg) or more or 8 years and above.	1 set	M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)
M5072A CE	Infant/Child SMART Pads Cartridge for infants/ children under 8 years old or 55 lbs. (25kg).	1 set	M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)
989803139261 CE	SMART Pads II for adults, infants and children.	1 set	861304, 861305 (HeartStart FRx Defibrillators)



M3870A



M5071A



M5072A



989803139261

AED Batteries


Product No.	Description	Shelf Life	Philips Instrument
BT1 CE	ForeRunner lithium manganese dioxide battery.	5 year shelf life plus 1 year standby life	E, S, EM (Heartstream ForeRunner AED) 940010XX, 940020XX, 940030XX (Heartstart FR AED)
M3848A CE	FR2+ rechargeable battery	2.5 years shelf life 6 month standby life	M3860A, M3861A, M3840A, M3841A (HeartStart FR2+ AED sold after 12/01/01)
M3849A CE	Charger for the FR2+ rechargeable battery	na	
M3863A CE	FR2 series lithium manganese dioxide battery.	5 year shelf life plus typical 5 year standby life	M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 Series AED)
M5070A CE	Battery for HeartStart HS1/FRx defibrillators (lithium manganese dioxide).	5 year shelf life plus typical 4 year standby life	M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators) 861304, 861305 (HeartStart FRx Defibrillators)
989803136291 CE	Battery-Aviation-LiMnO2- FR2 series. Meets requirements of FAA TSO C-142 regulation for use on commercial aircraft	5 year shelf life plus typical 5 year standby life	M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 Series AED)
989803136301 CE	Battery-Aviation-LiMnO2-ForeRunner Meets requirements of FAA TSO C-142 regulation for use on commercial aircraft	5 year shelf life plus 1 year standby life	E, S, EM (Heartstream ForeRunner AED) 940010XX, 940020XX, 940030XX (Heartstart FR AED)



M3848A, M3849A



M3863A



M5070A

AED Accessories and Kits

Product No.	Description	Philips Instrument
M3873A AAMI	FR2+ ECG Assessment Module. 3-wire ECG module lets you use FR2+ with ECG display to assess ECG of responsive, breathing, but at-risk patient. (AAMI lead colors).	M3840A, M3860A (HeartStart FR2+ AED with ECG sold after 12/01/01)
M3874A CE IEC	FR2+ ECG Assessment Module. 3-wire ECG module lets you use FR2 series with ECG display to assess ECG of responsive, breathing, but at-risk patient. (IEC lead colors).	
68-PCHAT	AED Fast Response Kit. Kit includes two pairs of nitrile gloves, one CPR mask, one paramedic scissors, one prep razor, one large reinforced paper towel. All items contained in a zippered pouch that attaches to the AED carrying case.	E, S, EM (Heartstream ForeRunner AED) 940010XX, 940020XX, 940030XX (Heartstart FR AED) M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 series AED) M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators) 861304, 861305 (HeartStart FRx Defibrillator)
989803139311 CE	Infant/Child Key used in infant/child defibrillation with the HeartStart FRx Defibrillator. For infants and children under 8 years old or <55 lbs. (<25kg).	861304, 861305 (HeartStart FRx Defibrillator)

AED Electrode Pad Adapters

- for cross-product compatibility between FR/FR2/FRx series pads and manual defibrillators


Product No.	Description	Defibrillator Models
05-10000 CE	AED QUIK-COMBO electrode adapter	Medtronic Physio-Control LifePak 9, 10c, 11, 12, 500, 20
05-10100 CE	AED Zoll™ electrode adapter	Zoll 1200, 1400, 1600, and M-Series
05-10200 CE	AED Philips electrode adapter	Philips CodeMaster 100, XE, XL, XL+; Philips Heartstream/HeartStart XL and XLT; Laerdal Heartstart 4000 (Includes: Hewlett-Packard/Agilent branded versions)



M3873A



989803139311



05-10000, 05-10200, 05-10100
Pads Adapters

Data Management Solutions

- Helps improve the quality of your response to SCA through retrospective data analysis of previous responses

AED Data Cards and Accessories

Product No.	Description	Philips Instrument
M3524A	FR2 series data card reader for a personal computer.	M3834A (HeartStart Event Review Software) 989803141811 (HeartStart Event Review Organization-wide License) 861276 (HeartStart Event Review Pro Software) 861311 (HeartStart Review Express Connect Software) M3831A, M3833A (CodeRunner Web Software) M3854A (FR2 series data card and tray) M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 series AED)
M3853A	FR2 series AED tray for data card (M3854A). Houses FR2 series AED data card when inserted into AED.	
M3854A	FR2 series AED data card and tray. Stores device history and setup data.	M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 series AED)
SDCF05	Adapter, Data Card to PCMCIA. Enables the FR2 series data card (M3854A) to be read on a PC that has a PCMCIA format data card reader.	
ACT-IR	Infrared Data Cable for HeartStart AEDs. Enables data collection from HS1 or FRx to PC running Event Review Software.	M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators) 861304, 861305 (HeartStart FRx Defibrillators)

HeartStart Event Review Software


Product No.	Description	Philips Instrument
861276 CE	HeartStart Event Review Pro Software Full featured software that allows you to review and create reports on sudden cardiac arrest event data recorded by Philips defibrillators (English, French, German, Italian, Japanese, Spanish). A01 = Single PC License A02 = 3-PC License A03 = Site License	M3854A, M3524A (FR2 series data card and data card reader) M3500B (Heartstream/HeartStart XLT) M4735A (Heartstream/HeartStart XL) M3535A, M3536A, M3536J (HeartStart MRx, HeartStart MRxE)
M3834A CE	HeartStart Event Review Software Software that allows you to review and create reports on sudden cardiac arrest event data recorded by Philips defibrillators (English, French, German, Italian, Japanese, Spanish) for single PC license.	M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 series AED) E, S, EM (Heartstream ForeRunner AED)
989803141811	HeartStart Event Review Software Software that allows you to review and create reports on sudden cardiac arrest event data recorded by Philips defibrillators (English, French, German, Italian, Japanese, Spanish) for unlimited PC's in organization.	940010XX, 940020XX, 940030XX (Heartstart FR AED) 861304, 861305 (HeartStart FRx Defibrillators)
861311 A01	HeartStart Review Express Connect Software Allows you to download defibrillator data, view it, save it, and send to colleagues who have Event Review software via email. (English, French, German, Italian, Japanese, Spanish)	M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)
989803143041 CE	HeartStart Configure-PDA Software Allows you to configure your FRx or HS1 defibrillator (for use by personnel authorized by your medical director). Palm™ and Pocket PC. (English, French, German, Italian, Japanese, Spanish)	M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)
989803143051	HeartStart Case Capture-PDA Software Allows you to download case data from FRx or HS1 without removing AED from service. Palm™ OS. (English, French, German, Italian, Spanish)	861304, 861305 (HeartStart FRx Defibrillator)

AED Cases and Wall Mount Solutions


- Lets you place AEDs strategically and transport them to the patient quickly for fast response to a SCA emergency.

AED Wall Mounts/Accessories

Product No.	Description		Philips Instrument
M3857A	Wall bracket	Wall-mount the AED where it can be deployed quickly during a cardiac emergency.	E, S, EM (Heartstream ForeRunner AED) 940010XX, 940020XX, 940030XX (Heartstart FR AED)
M3858A	AED wall sign	Attaches to wall and extends out. Gives high visibility indication of AED location.	M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 series AED)
M3859A	Secure pull seal for wall mount (10 pack)	Plastic tie used to secure AED to wall mount bracket. Seal deters theft and indicates tampering.	861304, 861305 (HeartStart FRx Defibrillator) M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)



M3857A



M3858A



M3859A

AED Cases/Wall Cabinets

Product No.	Description	Dimensions	Philips Instrument
M3868A	Carrying case (red) for FR2 series AED (corduroy fabric).	Dimensions: 11" (W) x 9" (H) x 5" (D) 1 lb.	M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 series AED)
M3869A	Carrying case (red) for FR2 series AED (vinyl).		
989803136531	Defibrillator Cabinet-Basic. Basic cabinet with audible alarm. Accomodates ForeRunner, FR2 series, HS1, FRx as well as Fast Response Kit.	16.5" (W) x 15" (H) x 6" (D) 11 lbs.	M3860A, M3861A, M3840A, M3841A, (Heartstream/HeartStart FR2 series AED) E,S,EM (Heartstream ForeRunner AED)
PFE7023D	AED Cabinet, semi-recessed, with audible alarm and light. Can connect to central security systems.	Inside Dimensions: 14" (W) x 22.5" (H) x 6" (D) 21 lbs.	940010XX, 940020XX, 940030XX (Heartstart FR AED) M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)
PFE7024D	AED Cabinet, wall surface mounted, with audible alarm and light. Can connect to central security sytstems.	Dimensions: 16" (W) x 22.5" (H) x 6" (D) 20 lbs.	861304, 861305 (HeartStart FRx Defibrillators)
YC	Carrying case (red), hard plastic waterproof shell.	Dimensions: 12 1/4" (W) x 11 5/8" (H) x 5 3/4" (D) 4 lbs.	M3860A, M3861A, M3840A, M3841A, (Heartstream/HeartStart FR2 series AED) 861304, 861305 HeartStart FRx Defibrillators) M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)
M5075A	Standard carrying case for HeartStart HS1.	9.5" (W) x 8.5" (H) x 4.8" (D) 1.2 lbs.	M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)
M5076A	Slim carrying case for HeartStart HS1.	9.5" (W) x 8.5" (H) x 3.5" (D) 1 lb.	
989803133171 CE	Temperature Control Carrying Case. Maintains AED at standby temperature in vehicles in extreme cold and hot climates.	15 7/16" (W) x 7 7/16" (H) x 5" (D) 11.2 lbs.	M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 series AEDs)
989803139251	FRx carrying case made of wipeable urethane.	9.5" (W) x 8" (H) x 5" (D)	861304, 861305 (HeartStart FRx Defibrillators)



M3868A



M5076A




989803139251


AED Training Materials

- Helps your responders practice delivering care to SCA victims correctly and confidently

Product No.	Description	Quantity	Philips Instrument
07-10900 CE	Reusable AED adult training electrode pads with yellow connector. Non-clinical training and demonstration pads.	1 set	E, S, EM (Heartstream ForeRunner AED when equipped with Training Card) 940010XX, 940020XX, 940030XX (Heartstart FR AED equipped with Training Card) M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 series AED equipped with M3864A Training and Administration Pack) M3752A (AED Trainer 2) 07-10801 (AED Trainer)
07-11000	AED Trainer 2 carrying case	1	07-10801 (AED Trainer) M3752A (AED Trainer 2)
M3752A CE	AED Trainer 2. AED trainer that looks and behaves like FR2 series.	1	M3753A, M3754A, 07-10900, M3871A, M3755A, M3756A, 07-11000
M3753A	Remote control for AED Trainer 2	1	M3752A, M3756A (AED Trainer 2, AED Little Anne Training System)
M3754A	Optional programming kit for AED Trainer 2 (includes PC cable with software)	1	M3752A, M3756A (AED Trainer 2, AED Little Anne Training System)
M3755A CE	Training Pads for AED Little Anne. Works with AED Trainer 2 and AED Little Anne manikin to give the student realistic feedback on pads placement.	1	M3756A (AED Little Anne Training System)
M3756A CE	AED Little Anne Training System. All-in-one FR2 Series training solution. Contains AED Trainer 2, remote control, special training pads, and a manikin.	1	M3755A (Training Pads for AED Little Anne)
M3855A CE	Charger for FR2 series AED training and administration pack	1	M3864A (Training and Administration Pack for FR2 Series AEDs)
M3871A CE	HeartStart FR2 series Infant/Child Training Pads	1 set	M3864A, M3752A (Training and Administration Pack for FR2 Series AEDs, AED Trainer 2)



07-10900



M3871A

AED Training Materials (Continued)


Product No.	Description	Quantity	Philips Instrument
M3860-97800	Quick Reference Card, HeartStart FR2+ (English) Other languages available. Contact your order administrator.	1	M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 series AEDs)
M3864A CE	Training and Administration Pack for HeartStart FR2 AED. Allows AED to be configured and used for training without depleting AED battery. Rechargeable.	1	M3860A, M3861A, M3840A, M3841A, (Heartstream/HeartStart FR2 series AEDs) M3855A (Charger for Training and Administration Pack)
M5073A CE	HS1 Adult Training Pads Cartridge and Pad Placement Guide.	1	M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)
M5074A CE	HS1 Infant/Child Training Pads Cartridge and Pad Placement Guide.	1	
M5093A CE	HS1 Replacement Adult Training Pads with replaceable seal.	1	M5073A (Adult Training Pads Cartridge and Pad Placement Guide)
M5094A CE	HS1 Replacement Infant/Child Training Pads with replaceable seal.	1	M5074A (Infant/Child Training Pads Cartridge and Pad Placement Guide)
M5066-89100 CE	HS1 Instructor's Training Toolkit NTSC (PAL is M5066-89101).	1	M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)
989803150291	FR2+ Instructor's Training Toolkit NTSC (PAL is 989803150301)	1	M3860A, M3861A, M3840A, M3841A (Heartstream/HeartStart FR2 series AEDs)
M5066-97800 CE	Quick Reference for HeartStart HS1 (English) Other languages available. Contact your order administrator.	1	M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)
989803139271 CE	Training Pads II facilitates training on the HeartStart FRx Defibrillator. Can be used to train for adult and infant/child resuscitation (when used with Infant/Child Key).	1	861304, 861305 (HeartStart FRx Defibrillators)



M3864A



M5073A



M5074A



989803139271

AED Training Materials (Continued)

Product No.	Description	Quantity	Philips AED Instrument
M5085A CE	HeartStart Trainer. Resembles HeartStart HS1 Defibrillator. Provides eight realistic training scenarios and comes with nylon carrying case, one set of Adult Training Pads and one External Manikin Adapter (English).	1	M5066A, M5067A, M5068A (HeartStart HS1 Defibrillators)
M5088A	Internal Manikin Adapter for use with HeartStart Trainer. Compatible with Laerdal ResusciAnne and Little Anne CPR manikins.	1	M5085A (HeartStart Trainer) 861306 (HeartStart FRx Trainer)
M5089A	External Manikin Adapter for use with HeartStart Trainer and all CPR manikins.	5/roll	M5085A (HeartStart Trainer) 861306 (HeartStart FRx Trainer) All CPR manikins
861306	HeartStart FRx Trainer. Resembles HeartStart FRx Defibrillator. Provides eight realistic training scenarios and comes with nylon carrying case, one set of Training Pads II and one External Manikin Adapter. (U.S. English, UK English, French, German)	1	861304, 861305 (HeartStart FRx Defibrillators)



M5085A



861306

Numerics

05-10000 Pad adapter (Medtronic Physio Control), 3, , 29
05-10100 Pad adapter (Zoll), 3, , 29
05-10200 adapter (HP, Agilent, Philips, Laerdal), 3, , 29
07-10900 Reusable AED adult training electrode pads, 34
07-11000 AED Trainer 2 carrying case, 34
40401A Reusable Infant NIBP Cuffs, 16
40401B Reusable Pediatric NIBP Cuffs, 16
40401C Reusable Adult NIBP Cuffs, 16
40401D Reusable Large Adult NIBP Cuffs, 16
40401E Reusable Thigh NIBP Cuffs, 16
40457C Chemical/thermal recording paper, 11
40457D Chemical/thermal recording paper, 11
68-PCHAT AED Fast Response Kit, 29
861276 HeartStart Event Review Pro, 13
861276 HS Event Review Pro Software, 31
861306 HeartStart FRx Trainer, 36
861311 A01 HeartStart Review Express Software, 31
9898031 47711 HeartStart data card, 13
989803133171 Temp. Control Case, 33
989803136291 Battery-Aviation-LiMnO₂, 28
989803136301 Battery-Aviation-LiMnO₂, 28
989803136531 Defib Cabinet-Basic, 33
989803138171 Recording paper, 11
989803138181 Recording paper, 11
989803139251 FRx Case, 33
989803139261 SMART Pads II, 27
989803139271 Training Pads II, 35
989803139311 Infant/Child Key, 29
989803141811 HS Event Review Software, 31
989803142521 HeartStart transfer station, 13
989803143041 HeartStart PDA Software, 31
989803143051 HeartStart Capture-PDA Software, 31
989803150291 FR2+ Instructor's Training Toolkit NTSC, 35

A

ACT-IR, Infrared data cable, HS1, 30
Adapter Cable Test Loads, 3

B

BT1 ForeRunner lithium manganese dioxide battery, 28

C

Continuous Temperature

21075A Reusable Esophageal/rectal probe, 25
21076A Reusable Esophageal/rectal probe, 25
21078A Reusable skin surface probe, 25
21082A Extension cable for disposable probes, 24
21082B Extension cable for disposable probes, 24
21090A Disposable temperature probe 12FR, 24
21091A Disposable temperature probe, 24
21093A Disposable temperature probe 12FR, 24

21094A Disposable temperature probe 18FR, 24
21095A Disposable temperature probe 24FR, 24
21096A Disposable temperature probe 16FR, 24
21097A Disposable temperature probe 18FR, 24
M1837A Disposable temperature probe 9FR, 24
M2255A Disposable temperature probe 14FR, 24

D

Disposable adult electrode pad, 1
Disposable pediatric electrode pad, 1
Disposable pediatric radiotransparent electrode pad, 1
DP2 Self-adhesive adult electrode pads, 27
DP6 Self-adhesive adult electrode pads, 27

I

Invasive Pressure

CPJ84022 Disposable domes, 23
CPJ84046 Transducer holder, 23
CPJ840J6 Reusable pressure transducer, 23
M1567A Disposable BP kit, 23
M1568A Disposable BP kit, 23
M1634A adapter cable, 23

M

M1131A Adult/Pediatric disposable sensor, 20
M1191A SpO₂ Sensor, Adult Finger, 19
M1192A SpO₂ Sensor, Pediatric/Adult Finger, 19
M1194A Adult ear clip SpO₂ sensor, 19
M1195A SpO₂ Infant sensor, 19
M1572A Comfort NIBP cuff, Pediatric, 16
M1573A Comfort NIBP cuff, Small Adult, 16
M1574A Comfort NIBP cuff, Adult, 16
M1575A Comfort NIBP cuff, Large Adult, 16
M1598B NIBP Interconnect Tubing, 18
M1599B NIBP Interconnect Tubing, 18
M1603A 3-lead grabber, 7
M1605A 3-lead snap, 7
M1613A 3-lead grabber, 7
M1613A lead set, 7
M1615A 3-lead snap, 7
M1623A 5-lead grabber, 8
M1625A 5-lead snap, 8
M1633A 5-lead grabber, 8
M1635A 5-lead snap, 8
M1644A 5 lead snap, 8
M1645A 5 lead snap, 8
M1672A Lead set, 7
M1673A Lead set, 7
M1674A Lead set, 7
M1741A Switchless Internal paddle with one-piece handle and electrode, 5
M1742A Switchless Internal paddle with one-piece handle and electrode, 5

- M1743A Switchless Internal paddle with one-piece handle and electrode, 5
- M1744A Switchless Internal paddle with one-piece handle and electrode, 5
- M1781A CodeMaster test load, 3
- M1783A Sync cable, 6
- M1789A Replacement paddle electrode, 4
- M1920A Microstream FilterLine Adult/Ped, 21
- M1921A Microstream FilterLine H Adult/Ped, 21
- M1923A Microstream FilterLine H Inf/Neon., 21
- M1941A Extension Cable for Reusable SpO₂ Sensors, 19
- M1943A Sensor adapter cable, 20
- M1968A 5 lead grabber, 8
- M1971A 5 lead grabber, 8
- M2202A Disposable electrode, 11
- M2520A Microstream Smart CapnoLine O₂ Ped., 21
- M2522A Microstream Smart CapnoLine Plus O₂ Adult, 21
- M2524A Microstream Smart CapnoLine Ped., 21
- M2526A Microstream Smart CapnoLine Plus Adult, 21
- M3501A Philips (AAMI) adult electrode pads, 1
- M3502A Philips (IEC) adult electrode pads, 1
- M3503A Philips (IEC) pediatric electrode pads, 1
- M3504A Philips (AAMI) pediatric electrode pads, 1
- M3506A HeartStart battery charger adapter, 12
- M3507A Multi-function defibrillator electrode cable, 2
- M3508A HeartStart defibrillator hands-free electrode cable, 2
- M3508A HeartStart hands-free pad cable, 2
- M3509A HeartStart carrying case, 12
- M3516A HeartStart battery, 12
- M3524A FR2 data card reader, 30
- M3537A HeartStart MRx bed rail hook, 14
- M3538A HeartStart MRx battery, 12
- M3541A HeartStart MRx carrying case, 14
- M3542A External anterior/anterior paddle set, 4
- M3543A External Anterior/Anterior Paddle Set, 4
- M3544A HeartStart data card tray, 14
- M3545A HeartStart MRx data card and tray, 14
- M3713A HeartStart Adult Plus electrode pads, 1
- M3716A HeartStart Adult Radiolucent electrode pads, 1
- M3717A HeartStart Pediatric Plus electrode pads, 1
- M3718A HeartStart Adult Radiopaque/Reduced Skin Irritation electrode pads, 1
- M3719A HeartStart Pediatric Radiopaque/Reduced Skin Irritation electrode pad, 1
- M3725A HeartStart test load, 3
- M3752A AED Trainer 2, 34
- M3753A Remote control for AED Trainer 2, 34
- M3754A Optional programming kit for AED Trainer 2, 34
- M3755A Training Pads for AED Little Anne, 34
- M3756A AED Little Anne Training System, 34
- M3834A HS Event Review Software, 31
- M3848A FR2+ rechargeable battery, 28
- M3849A Charger for the FR2 rechargeable battery, 28
- M3853A FR2 AED tray for data card M3854A, 30
- M3854A FR2 AED data card and tray, 30
- M3855A Charger for FR2 AED training and admin pack, 34
- M3857A Wall mount for FR2 AED, 32
- M3858A AED wall sign, 32
- M3859A Secure pull seal for wall mount, 32
- M3860-97800 Quick Reference Card, HS FR2+, English, 35
- M3863A FR2 lithium manganese dioxide battery, 28
- M3864A Training and Admin pack for HS FR2 AED, 35
- M3868A Carrying case for Heartstream FR2 AED, 33
- M3869A Carrying case for FR2 AED, 33
- M3870A Self-adhesive FR2 electrode pads, 27
- M3871A HeartStart FR2 Training pads, 34
- M3873A FR2+ ECG Assessment Module, 29
- M3874A FR2+ ECG Assessment Module, 29
- M4552A Infant, Antimicrobial Coated NIBP Cuff, 15
- M4553A Pediatric, Antimicrobial Coated NIBP Cuff, 15
- M4554A Small Adult, Antimicrobial Coated NIBP Cuff, 15
- M4555A Adult, Antimicrobial Coated NIBP Cuff, 15
- M4557A Large Adult, Antimicrobial Coated NIBP Cuff, 15
- M4559A Thigh, Antimicrobial Coated NIBP Cuff, 15
- M4572A Infant, Disposable NIBP Cuff, 17
- M4573A Pediatric, Disposable NIBP Cuff, 17
- M4574A Small Adult, Disposable NIBP Cuff, 17
- M4575A Adult, Disposable NIBP Cuff, 17
- M4576A Extra Long Adult, Disposable NIBP Cuff, 17
- M4577A Large Adult, Disposable NIBP Cuff, 17
- M4578A Extra Long, Large Adult, Disposable NIBP Cuff, 17
- M4579A Thigh, Disposable NIBP Cuff, 17
- M4740A HeartStart Switchless Internal paddle adapter, 6
- M4741A Switched Internal paddle electrodes, 5
- M4742A Switched Internal paddle electrodes, 5
- M4743A Switched Internal paddle electrodes, 5
- M4744A Switched Internal paddle electrodes, 5
- M4745A External anterior/anterior paddle set, 4
- M4746A External anterior/anterior paddle set, 4
- M4748A Adapter extension cable for Heartstream XL Defibrillator/Monitor, 6
- M4748A HeartStart extension cable, 2
- M4751A Accessory pouch for Heartstart XL, 12
- M4759A Paddle Electrode Replacement, 4
- M5066-89100 FR2+ & HS1 Instructors Training Toolkit, 35
- M5066-97800 Quick Reference for HS1, 35
- M5070A HS1/FRx battery, 28
- M5071A Adult SMART Pads Cartridge, 27
- M5072A Infant/Child SMART Pads Cartridge, 27
- M5073A HS1 Adult Training Pads Cartridge, 35
- M5074A HS1 Infant/Child Training Pads Cartridge, 35
- M5075A Slim carrying case for HS1, 33
- M5076A Slim carry case for HS1, 33
- M5085A HeartStart Trainer, 36
- M5088A Internal Manikin Adapter, 36
- M5089A External Manikin Adapter, 36
- M5093A HS1 Replacement Adult Training Pads, 35
- M5094A HS1 Replacement Infant/Child Training Pads, 35
- M5509A Carrying case, 12
- M5521A-M5525A HS MRx colored handles, 14

M5526A Sync Cable, 6
 M5528A HeartStart MRx mounting bracket, 14
 M5529A HeartStart MRx DC Power Module, 14
 Multifunction defibrillator electrode adapter cable, 2

P

Patient Cables and Leads

10-Electrode Cable Sets, 9

M1663A 10-Electrode trunk cable, 9
 M1968A ICU grabber limb, 9
 M1971A ICU grabber limb, 9
 M1976A ICU grabber chest, 9
 M1978A ICU grabber chest, 9
 M3525A 10-Electrode trunk cable, 10
 M3526A 3-wire snap leadset, 10
 M3527A 7-wire snap leadset, 10
 M3528A 3-wire snap leadset, 10
 M3529A 7-wire snap leadset, 10
 M5530A Combiner plug, 10

3-Lead Monitoring, 7

M1500A Trunk cable leadset, 7
 M1510A Trunk cable leadset, 7
 M1580A Trunk cable leadset, 7
 M1590A Trunk cable leadset, 7
 M1669A Trunk cable leadset, 7
 M1733A Trunk cable leadset, 7
 M1735A Trunk cable leadset, 7

5-Lead Monitoring, 8

M1520A Trunk cable, 8
 M1530A Trunk cable, 8
 M1600A Trunk cable leadset, 8
 M1610A Trunk cable leadset, 8
 M1668A Trunk cable leadset, 8
 M1734A Trunk cable leadset, 8
 M1736A Trunk cable leadset, 8

PFE7023D AED Cabinet, semi-recessed, with audible alarm and light, 33

PFE7024D AED Cabinet, wall surface mounted, with audible alarm and light, 33

R

Recording Paper

Chemical/Thermal Recording Paper, 11

S

SDCF05 Adapter, Data Card to PCMCIA, 30

Y

YC Carrying case, hard plastic waterproof shell, 33



Philips Medical Supplies

Philips is committed to producing and supporting the finest quality medical equipment and supplies. Our supplies are thoughtfully designed, tested and manufactured to deliver reliable and accurate results from your HeartStart defibrillator/monitors and AED's. For a complete list of supplies, please visit <http://shop.medical.philips.com>.

Philips Medical Systems is part of
Royal Philips Electronics

On the web:

<http://philips.com/heartstart>

Via e-mail:

medical@philips.com

By fax:

+31 40 27 64 887

By Postal service:

Philips Medical Systems
3000 Minuteman Road
Andover, MA 01810-1085

Asia:

Tel: +852 2821 5888

Europe, Middle East, Africa:

Tel: +31 40 27 87246

Latin America:

Tel: +1 954 628 1000

North America:

Tel: +1 800 285 5585

© Koninklijke Philips Electronics N.V. 2007

All rights are reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder.

Philips Medical Systems Netherland B.V. reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in the Netherlands

4522 962 21671/864 JAN 2007